Operating instructions

Powerful compact device for the determination of total organic carbon (TOC) in the trace range

The information in this manual is subject to change without notice and should not be construed as a commitment by the manufacturer. The membraPure GmbH accepts no responsibility for any errors that may appear in this manual.

Unless expressly stated otherwise in this manual, nothing in this manual shall be construed as a guarantee or warranty against loss, damage to persons or property, fitness for a particular purpose or the like.

The membraPure GmbH shall in no event be liable for incidental or consequential damages arising out of the use of the contents of this manual or the products described herein.

This manual or parts thereof may not be used without the written permission of membraPure GmbH may not be reproduced or copied and its content may not be passed on to third parties or used for unauthorized purposes. Violations will be prosecuted. Additional copies of this manual can be requested from membraPure GmbH at the applicable price.

© **Copyright 2022-2025** membraPure GmbH
All rights reserved.

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Contents

I	English	7
1	Product overview	9
1.1	Product overview	10
1.2	System architecture	10
2	Device and documentation instructions	11
3	Safety instructions	13
3.1	Occupational safety	14
3.2	Operational safety	14
3.3	Personnel requirements	14
4	Product installation	15
4.1	Goods receipt and inspection	16
4.2	Installation in a real environment	17
4.3	Installation of the glas reactor	18
4.4	Installing the UV lamp	21
4.5	Power connection	24
4.6	Connection options for samples and devices	24
4.6.1	Connectivity options on the device	25
4.6.2	Analog inputs and outputs	26
4.6.3	Autosampler	27
5	Operating options	29
6	Configuration	31
7	GUI - Graphical User Interface	33
8	Commissioning the device	37
8.1	Setting the time zone	38
8.2	Rights management	39
8.3	Create user	39
8.4	CFR option in user management	43
9	Operation	45
9.1	Login	46
9.2	Monitor	46
9.2.1	Warnings	49
9.2.2	UV intensity too low	49
9.2.3	Leakage in the system	49
9.2.4	Multiline-Sequence	49
9.3	Sample	51
9.3.1	Save Template	53
9.3.2	Load Template	54
9.4	Qualification	56
9.4.1	Calibration TOC	56
9.4.2	Calibration-Sequence	58
9.4.3	SST	63
9.4.4	Calibration of conductivity sensors cell constant	65
9.5	Direct	68
9.6	About	71
9.7	Status	72
9.8	Measurement data backup and visualization	73
9.8.1	Cloud-based export of data	73
9.8.2	Export measurement data	74
9.8.3	Display and report output of data	75
9.9	DataViewer	76
9.9.1	Data input dialogue	77
9.9.2	Date selection dialogue	80
9.9.3	Data display dialogue	81
9.9.4	Document preview	85
9.9.5	Data export dialogue	86
9.10	Remote access	87
9.11	Messages	88
9.12	Shut down system	89
10	Settings	91
10.1	Quick setup	92
10.2	Device settings	96
10.2.1	Change UV	96
10.2.2	Change Pump Tube	96
10.2.3	Backup	97
10.2.4	Update Software	98
10.2.5	License	100
11	Factory settings	101
11.1	User reset	102
11.2	Rights reset	103
11.3	Settings reset	103
11.4	Device factory settings	104
12	Repair and service	107
12.1	Cloud	108

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Product overview

1.1. Product overview

The system is used to analyze the TOC content in aqueous solutions. The organic substances are oxidized using short-wave and high-energy UV radiation. The resulting CO_2 forms carbonic acid and hydrogen carbonate on contact with an aqueous solution. The resulting charges in the solution lead to a measurable increase in conductivity. This conductivity is then measured with a pair of electrodes and converted into the TOC value.

1.2. System architecture


Figure 1.2.1.: System sequence with component description

This system has four (4) inputs (input: 2 to 5) for sample intake and one sample output (output: 1). The sample inputs can be controlled individually via three different solenoid valves. The user is free to control these inputs manually or have them controlled automatically by the system controller. The input conductivity of the sample is checked at LF1 (conductivity sensor 1) before oxidation. After the UV reactor, the conductivity of the oxidized sample is measured again and the TOC value of the sample can be determined via a previous calibration.

Chapter 2

Device and documentation instructions

Table 2.0.1.: Description of the notes in the document

Note	Meaning
	This notice draws your attention to situations that can lead to material damage.
	These instructions serve as a guide for operating the appliance and should, but do not have to be followed.
	This notice warns you that failure to follow the instructions may result in injury.
	This notice draws your attention to the fact that failure to follow the instructions may result in death or serious injury.

Table 2.0.2.: Description of the instructions on the device

Note	Meaning
	Risk of injury! UV radiation. Strongly oxidizing. May cause injuries to the skin and eyes!
	Risk of injury! Formation of ozone, can be harmful to cells.
	Danger to life! Warning of electric current!

Chapter 3

Safety instructions

3.1. Occupational safety

The user is responsible for complying with the following instructions:

- Installation instructions
- Local standards and instructions

3.2. Operational safety

Before using the appliance, it should be checked for damage. This includes

- External damage to the housing
- Damage to the sample lines
- Damage to the power cable or power supply unit and mains plug
- Damage to sample lines, these must not be kinked
- Sample lines must be hand-tightened

If damage is detected:

1. Do not put the device into operation
2. Inform the technical support of the sales service provider or the manufacturer directly
3. Sample lines that are too loose can be screwed on hand-tight

3.3. Personnel requirements

- Installation, commissioning, operation and maintenance may only be carried out by trained personnel.
- Training is carried out by means of e-learning or service training.
- Personnel must have read and understood the operating instructions.
- The instructions in the operating instructions must be followed.
- Faults on the appliance may only be rectified by trained personnel.

Faults that are not described may only be rectified directly by the manufacturer or by the service organization. For further information, please read chapter 3.2.

Chapter 4

Product installation

4.1. Goods receipt and inspection

- When accepting the delivery, pay attention to:
 - Damage to the packaging
 - Damage to the delivery
 - Completeness of the delivery
- The delivery includes separately supplied glass reactors and a quartz-mercury UVC lamp.
- If damage is detected:
 - Let the supplier know
 - Keep packaging until the case has been clarified
 - Report the damage to membraPure GmbH

4.2. Installation in a real environment


Figure 4.2.1.: Installation on a flat surface

Before installation, chapter 3 safety instructions must be read and understood! Otherwise, incorrect installation cannot be ruled out.

- Installation only on a firm, level surface
 - Load capacity: min. 30 kg
 - Surfaces free of dirt, oil and grease
- The product may only be installed vertically!
 - Front and side areas should be freely accessible
 - Minimum distance to other devices: min. 50 cm
 - Recommended distance: 1 m
- A seamless connection to the power grid is mandatory; refer to chapter 4.5
- When connecting to a water ring line system, the specified dimensions and connection interfaces must be strictly adhered to.
- Environmental requirements:
 - For indoor use only
 - Installation height: max. 3500 m
 - Ambient temperature: 10°C - 45°C
 - Relative humidity: 10 % - 95 %
- For the sample connection, please read chapter 4.6.

4.3. Installation of the glass reactor

Do not start the appliance without the UV protection cover installed! Operating the device without the UV protection cover can damage the device and, in the worst case, lead to injuries!

The glass reactor may only be installed with the utmost care and by trained personnel - risk of breakage!

For correct and safe installation of the glass reactor, please follow the instructions below. The installation requires the following parts:

- Glass reactor (externally packaged)
- Security key
- 3 mm Allen key (hexagon socket) - spanner

Installation steps:

1. Open the front door of the appliance with the enclosed security key


Figure 4.3.1.: Locking the front door

2. Open the upper reactor bracket by loosening the two 3 mm hexagon socket screws. Remove the front part of the retaining bracket of the reactor (1) and the retaining brackets of the tubes (2,3).


Figure 4.3.2.: loosen and open the upper reactor bracket (1) and tube brackets (2,3)

3. Place the enclosed glass spiral reactor in the lower reactor holder (5). **Attention: Fragile!**


Figure 4.3.3.: Insert the reactor into the holder (5) and attach the upper reactor holder (4)

4. Attach the glass spiral reactor to the holder (4) - **the glass tube should still be able to move in the holder, this is not a mistake!!!**
5. Now connect the reactor connectors to metal fittings of the sample tubes. Fix the connections by hand (8).


Figure 4.3.4.: Screw the fittings (8) to the reactor glass coil

8. Now place the brackets of the two hose guides and tighten the screws (9).
9. During the first installation, you can now also install the UV lamp. To do this, go to the chapter 4.4 .
10. Install the UV protection cover with the appropriate screws. Ensure that the sample tubes pass through the cut-out on the right-hand side of the UV cover.

4.4. Installing the UV lamp

Do not start the device without the UV protective cover installed! Operation without the UV protective cover can damage the device and, in the worst case, lead to injuries!

The UV lamp may only be installed with the utmost care and by trained personnel - risk of breakage! Do not touch the UV lamp with bare hands.

The UV lamp may only be installed with the utmost care and by trained personnel - risk of breakage! Do not touch the UV lamp with bare hands. For correct and safe installation of the UV lamp, please follow the instructions below. No tools are required.

1. Open the front door of the appliance with the enclosed key.


Figure 4.4.1.: Locking the front door

2. Remove the installed UV protection cover by loosening the retaining screws (1)


Figure 4.4.2.: Screw connections of the UV cover

3. Loosen the retaining screw of the fluidic plate (3). Fold down the entire fluidic plate slightly to insert the UV lamp into the lower holder.


Figure 4.4.3.: Folding down the fluidic plate

4. Place the UV lamp (6) from above into the installed glass spiral reactor. If this is not yet installed, please first read chapter 4.3. When inserting the UV lamp, the lower connecting cables and the plug connector must first be threaded through the glass filament reactor. **Attention: Fragile! Handle with care!**


Figure 4.4.4.: Inserting and connecting the UV lamp

5. Now insert the two plugs into the sockets provided (7).
6. Return the front panel of the appliance to its original position and secure it with the corresponding screws
7. Install the UV protective cover with the corresponding screws. Ensure that the sample tubes pass through the recess on the left-hand side and that the cables of the UV lamp are in the recesses on the right-hand side.

4.5. Power connection


Device energized while in operation. Danger to life due to electric shock! When the main switch is switched off, the mains filter, overvoltage module and main switch are live. Disconnect the appliance from the mains (pull out the mains plug).

- Before connecting, ensure that the mains voltage corresponds to the specified mains voltage.
- Before connecting, check that the mains connection complies with local and national regulations.
- The system requires a power socket.
- Safe operation can only be guaranteed if the following standard voltage is used:
 - Europe/Asia: 230V $\sim \pm 10\%$, 50 - 60 Hz
 - Australia: 240V $\sim$, 50 Hz
 - USA/Canada: 120V $\sim$, 60 Hz
- The system is supplied with European plugs. If the sockets in the respective country do not allow the connection of these European plugs, adapters are required.
- The model with internal power supply unit is earthed via the power supply unit. Additional earthing via the housing is possible.

4.6. Connection options for samples and devices

The electrical connections for devices and the connections for samples are located on the right-hand side of the device. The labeling on the device clearly identifies the electrical and fluid connections. Refer to the table 4.6.1 for details.

Number	Connection	Standard	Function
1	Power ON/OFF	-	Switching the device on/off
2	USB	2.0	Connection
3	Connector 1	6-pin - Analog OUT	Connection
4	Connector 2	6-pin - Analog OUT	Connection
5	LAN port	Ethernet connection	Connection
6	Autosampler	4-pin - Communication Interface	Connection
7	Port 1	UNF ¼ - 28"	Sample receipt
8	Port 2	UNF ¼ - 28"	Sample receipt
9	Port 3	UNF ¼ - 28"	Sample receipt
10	Port 4	UNF ¼ - 28"	Sample receipt
11	Waste	UNF ¼ - 28"	Sample output
12	Electricity	IEC power connection	Electricity

4.6.2. Analog inputs and outputs

The update interval for the analog outputs is one minute. During a calibration or the SST, the last value is displayed until a new measurement is started. The standby signal corresponds to NAMUR NE 43 and is selectable between 0 and 4 mA.

Assignment of the sockets


Figure 4.6.2.: Connector pin assignment

Table 4.6.2.: Pin assignment connector 1

Pin	Description
1: 4(0) - 20 mA (GND)	TOC (0 to set limit value / limit)
2: 4(0) - 20 mA (+)	
3: 4(0) - mA (GND)	Input conductivity (0 to set limit value / limit)
4: 4(0) - mA (+)	
5: Relay; potential-free	Collective fault message or limit value exceeded depending on the selected option
6: Relay; potential-free	

- Pin 1/2: supplies an analog signal between 4 and 20 mA for the TOC measured value at the Port 1 input between 0 and the limit value/limit (set under DIRECT MODE, IO 1)
- Pin 3/4: supplies an analogue signal between 4 and 20 mA for the input conductivity between 0 and the limit value/limit in $\mu\text{S}/\text{cm}$ (set under DIRECT MODE, IO 2)
- Pin 5/6: Normally open - system error and limit value violation depending on the selected option

The measured value is updated every minute. During a calibration or the SST, the last value is displayed until a new measurement is started. The standby signal corresponds to NAMUR NE 43 and can be selected between 0 and 4 mA.

Table 4.6.4.: Pin assignment connector 2

Pin	Description
1: 4(0) - 20 mA (GND)	TOC 2 (0 to set limit value / limit)
2: 4(0) - 20 mA (+)	
3: 4(0) - 20 mA (GND)	TOC 3 (0 to set limit value / limit)
4: 4(0) - 20 mA (+)	
5:Control (GND)	Control input / trigger (Relay, potential free)
6:Control (+)	

Pin 1/2: supplies an analog signal between 4 and 20 mA for the TOC measured value at input INPUT 2 between 0 and the limit value/limit (set under DIRECT MODE, IO 3)

Pin 3/4: supplies an analog signal between 4 and 20 mA for the TOC measured value at input INPUT 3 between 0 and the limit value/limit (set under DIRECT MODE, IO 4)

Pin 5/6: external control input / trigger, the measurement runs on closed relays .

The measured value is updated on-line. During a calibration or the SST, the last value is displayed until a new measurement is started. The standby signal corresponds to NAMUR NE 43 and can be selected between 0 and 4 mA.

4.6.3. Autosampler

The system offers the option of connecting an autosampler and operating the device in sample mode, 9.3. The scope of delivery of the sampler (option) includes:

- External power supply unit for the autosampler
- Serial connection cable (RS-232, 4-pin autosampler connection to 9-pin Sub-D socket).


Figure 4.6.3.: Arrangement of the connections on the rear panel of the autosampler

Connect the power supply unit to the autosampler (1). Ensure that the power supply is correctly connected and switched on. The autosampler's operating indicator (LED) can be used to check whether it is active.

Connect the serial connection cable:

- Plug and screw the 4-pin autosampler connection into the device
- Connect the serial connection cable (RS-232, 9-pin Sub-D socket) to the COM-1 port of the sampler (2).

Chapter 5

Operating options


Figure 5.0.1.: Device view

- The device is mainly operated via the 8" touchscreen (1)
- It is possible to connect external control elements.
 - A USB 2.0 port (2) is available on the right-hand side for this purpose (see chapter 4.6).
 - Keyboards, mice, USB sticks and printers can be connected there.
 - It is recommended to connect a USB HUB to be able to operate several devices simultaneously.
- The device software can be accessed via the screen
- The graphical user interface of the software is described in chapter 7 and chapter 9 .

Chapter 6

Configuration

- The device is available in the following versions:
- miniTOC Basic (inputs: sample and 1 inlet for calibration / SST)
 - miniTOC Pro (Inputs: Sample and 3 inlets for calibration / SST)
- The variants differ in the points:
- The Pro (Professional) version has 3 inputs for qualification and calibration. The Basic version has only one additional input.
 - The miniTOC Pro (Professional) is able to switch between port 1, port 3 and port 4 in multiline mode and thus measure different inputs sequentially.
- Depending on the options and variants selected, the appliances are equipped with the following components and features:
- Inlet pressure monitoring with alarm switch-off in the event of overpressure or missing supply.
 - Leakage monitoring with alarm function and locking of ports 1;3 and 4 at the input
 - Monitoring the system temperature
 - Autosampler and sample measurement mode for measuring individual samples
 - Extended analysis and logging functions for measurement and log data
 - Input and output connections for the adapted integration of the device into process systems

Chapter 7

GUI - Graphical User Interface

The meter has the following menus:

Table 7.0.1.: Description of the GUI

Monitor	Online measurement mode
Sample	Sample mode
Qualification	Calibration, SST (system suitability test), Conductivity
Data Viewer	Data output, Analysis, Reporting
Login	Login for users
Quick setup	Basic Settings
More	Advanced settings
Direct	Service mode
Remote	Remote connection
Service	Service
Status	Device status
About	About, Device Information
Help	Help manual


Figure 7.0.1.: Graphical user interface

Further menus are available via the More button on an additional screen page:

Table 7.0.2.: Description of the extended screen interface (More)

Settings	Device settings
Usermgm	User administration
Rights	Rights management
Multilines Sequence	Configuration of the multiline mode / definition of the sequence
Change UV	Confirm replacement of the pump hose of the UV-lamp
Change Pump Tube1	Confirm replacement of the pump hose of the pump1
Change Pump Tube2	Confirm replacement of the pump hose of the pump2
Export	Export of measurement data to an external data carrier
Update	Option to update the software
License	Option to extend the functional scope of the device by entering a license key
Change Timezone	Region-specific setting of the time zone
Shutdown	Shuts down the system to switch it off
Restart	Restarts the system including the operating system
Factory Reset	Restores the factory settings of the device settings, user settings or rights settings
One Time Service	Performing a service interaction, requires prior entry of a valid key
Backup	Provides various backup options: Measurement data, TOC software, system including operating system
Cloud	Provides the option of saving measurement data in the cloud via a network connection


Figure 7.0.2.: Advanced Graphical User Interface (More)

Chapter 8

Commissioning the device

8.1. Setting the time zone

To change the time zone, click on the button **More** and then to **Change Timezone**. You can now select your time zone.

If there is a time changeover between summer and winter time in your region, (Daylight Saving Time), activate this by clicking the checkbox.

Confirm the information with **OK** or reject the information with **Cancel**.

After changing the time zone, the software restarts automatically. The time adjusted to the new time zone is adopted after the restart and displayed in the main menu.


Figure 8.1.1.: Time zone pop-up menu

8.2. Rights management

The device offers you the option of dividing users into groups. Certain rights can be granted or withdrawn from these groups.

To set the user rights, click on the button in the main menu **More** . Then click on the button **Rights** .


Figure 8.2.1.: Rights management menu

8.3. Create user

After starting the device for the first time, you must create a user. Initially, only the service user and an administrator are created. These accesses are only available to authorized service specialists from the company membraPure GmbH and a user of rights group A (administrators).

The user of rights group A (administrators) is used to create new users on the system and is available after the first start of the system with the combination of user name and password:

- User name: Administrator
- Password: password

When logging in for the first time, the administrator's password must be changed, Figure 8.3.1.

Figure 8.3.1.: Change of password by the administrator on first login

To create a new user, click the button in the main menu **More**. Then please click on the button **Usermg**. To add a user, please click on the button **Add**.


The image shows the 'Create user menu' form. It has a title 'Please use only [0-9],[a-z],[A-Z]'. The form contains several fields: 'User-ID' (with a red border around 'New_User'), 'first name' (Hans), 'last name' (Mustermann), 'category' (group.B), 'is active' (checked), 'password' (masked with asterisks), 'pw limit' (UNLIMITED), 'active until' (X), 'confirm password' (masked with asterisks), 'new password at login' (checked), 'current_attempts' (0), 'unlimited attempts' (checkbox), 'Date created' (11.03.25), 'active from' (11.03.25), 'active until' (X), 'Active days' (90), 'remaining days' (X), and 'unlimited active' (checked). There are 'Add' and 'Cancel' buttons at the bottom, and 'Reset' buttons next to the 'unlimited attempts' and 'unlimited active' checkboxes.

Figure 8.3.2.: Create user menu

Please define a UserID, the group, first name and surname of the user. The assigned group also defines the user's authorizations. The table in Figure 8.3.3 shows all created rights groups.


Figure 8.3.3.: Rights groups are defined in the Rights management menu

To change the rights of the individual groups, please read chapter 8.2.

Also set the user's password. With the "New password at login" option, the user is prompted to choose a new password when logging in for the first time.

You can then see the user's current logins. With the "unlimited attempts" option, the password can be entered incorrectly any number of times.

Under "Active Days" you can specify how many days the user should exist. After the set time has expired, the account is irrevocably deleted. If you do not want to set a time limit, activate the "unlimited active" option on the right-hand side. If the entries for creating a new user are made using the **Add** the list of already created users appears in the result, Figure 8.3.4 .

After selecting the user in the list, the detailed information can be viewed using the **Edit** can be viewed or edited.

[illegible]

Figure 8.3.4.: Overview of existing users in the user administration

8.4. CFR option in user management

The CFR option is used to comply with regulatory requirements in accordance with 21 CFR Part 11 (Code of Federal Regulations), which are particularly applicable in regulated areas (e.g. pharmaceutical industry, laboratory environments).

- If the CFR option is activated, the software is configured so that it can uniquely identify users and log all actions in an audit-proof manner.
- If the CFR option is not activated, user administration is greatly simplified and does not meet the requirements of the above-mentioned regulation. This operating mode is not suitable for environments in which a clear assignment and traceability of user actions is required.

The CFR option can be purchased and activated at any time at a later date via a corresponding license.

If the CFR option has not been selected in the software, some restrictions apply. You can only log in by entering a user ID (e.g. user name or user ID). It is not necessary to enter a password. Users are not created or defined by the system. All users, with the exception of the administrative levels (administrator / service), receive the same authorizations. It is not possible to assign users to roles or authorization groups. No user groups can be defined.

In case the CFR option is activated, the behavior of the user administration changes. The software complies with the requirements for electronic records in accordance with 21 CFR Part 11. Each user must be created individually and assigned a unique user ID. A secure login with user name and password is mandatory. Users must be assigned to an authorization group. These groups control access to certain functions or data areas within the software. This means that user actions can be clearly assigned to a person and are logged to ensure complete traceability.

Chapter 9

Operation

9.1. Login

To log in as a user with your account, click on the button in the main menu **Login** . In the dialog box, (see Fig. 9.1.1) then enter your user name and password. Confirm the entry with **Ok** .

Figure 9.1.1.: Login pop-up menu

9.2. Monitor

Before using the device for the first time, it must be flushed for at least 30 minutes or until the same values are constantly displayed.

To start the measurement, connect the test sample to the device. Check that the UNF connection is hand-tightened.

Click on the button in the main menu **Monitor** . Log in with your user name and password.


Figure 9.2.1.: Monitor measurement mode selection menu

Now select the port that you want to measure online. With the **Singleline** -mode, only one of the four available ports can be measured at a time. To measure all ports at intervals, please select **Multiline**.


Figure 9.2.2.: Monitor measurement mode User interface

Number	Display	Description
1	Table field	Displays all measurements, including concentration [ppb], type, port, conductivity [$\mu\text{S}/\text{cm}$], time
2	Result field	Displays the result of the measurement, including factor, offset, coefficient of determination R2 and the percentage deviation from the current value
3	Info field	Displays information on the current measurement.
4	Login	Opens the login window
5	Comment	Opens the comment window
6	Report	Creates a report with the measurement data
7	Status	Opens the device status overview
8	Stop	Stops the measurement

You can use the info field to view useful information about the current measurement.

- Current Conc.(ppb) - Solution concentration of the current measurement sample
- Next Step - Next measurement step
- Status - Device status [Idle = idle, Flush = flush, measurement = measurement]
- Elapsed-Time - Displays the elapsed time
- User - Displays the logged in user
- ID - displays the measurement ID of the current measurement

After starting the measurement, the system is flushed for 300s. If you want to adjust the time, please read chapter 10.1.

If you want to cancel the measurement, press the Stop -button. Then enter your user name and password. You must then add a comment to the measurement (min. 4 characters). The measurement is then stopped and you are taken to the main menu.


Figure 9.2.3.: Pop-up menu Comment

9.2.1. Warnings

The device is able to display a warning as soon as the conductivity exceeds a threshold value. This threshold value must be set in the "Quick setup", as described in (chapter 10.1).

If the threshold value is exceeded, the TOC value and the conductivity are displayed in red digits and a warning appears in the info field.

9.2.2. UV intensity too low

This warning is displayed if the intensity of the UV lamp is too low or the UV lamp is defective. The device then automatically interrupts the measurement and a warning appears in the info field. If this occurs, the UV lamp must be replaced. Please contact our service staff (see chapter 12)

9.2.3. Leakage in the system

If there is a leak in the system, the device automatically stops the measurement and closes inlet valve 1. Once the leak has been eliminated, the measurement can be restarted. If the "Continue after Error" setting has been activated in the settings on the Options tab, the measurement is continued automatically.

9.2.4. Multiline-Sequence

The device offers the multiline sequence mode to measure several samples at certain intervals in the specified sequence. The intervals can be adjusted in the settings under

More and Multiline Sequence . In this way, the order of the ports and the associated measuring time can be set.

	Input		Time (min.)		
Step1	Port_1	▼	-	5	+
Step2	Port_3	▼	-	5	+
Step3	Port_4	▼	-	5	+
Step4	NA	▼	-	10	+
Step5	NA	▼	-	10	+
Step6	NA	▼	-	10	+
Step7	NA	▼	-	10	+
Step8	NA	▼	-	10	+

OK

Cancel

Figure 9.2.4.: Multiline sequence menu

Click on one of the drop-down menus under "Port". Then select one of the four available ports. You can then set the time to be measured under "Time (min.)".


Figure 9.2.5.: Multiline sequence menu Port selection

Confirm the entry with **OK**. You can then switch to monitor measurement mode. For more information on monitor measurement mode, please read chapter 9.2.

9.3. Sample

For easier operation, we recommend setting up the sample mode using a mouse and keyboard. For the correct connection of external input devices, please read chapter5.

Sample mode is used to measure samples in specific sequences. This sequence can be customized, saved and loaded.

To access the sample mode, navigate to the main menu and click on the button **Sample**. If you are not yet logged in, the device will ask you to log in with your user name and password. The graphical user interface will then open.


Figure 9.3.1.: Sample mode user interface

To set a measurement sequence, please follow the instructions below:

1. Select an item. To do this, click on "Position" in the corresponding line. A tab opens from which you can select various items


Figure 9.3.2.: Sample mode Position selection

2. Select a sample name. To do this, click on the "Designation" field and give the item a unique name (name, concentration). Example: Sucrose 500 ppb.


Figure 9.3.3.: Sample mode Sample name

3. "Date", "Time", "Cond" [Conductivity], "TOC" [TOC value] are filled in accordingly by the device during the measurement.
4. If you want to change the order of the measurements afterwards, you can do this using the buttons on the right-hand side.
5. If you want to delete an entry, click in the desired line and press the Delete button on the right-hand side.

9.3.1. Save Template

The device has the option to create and load templates for corresponding measurement sequences. To create a template, please follow the steps below:

1. Enter your measurement sequence and descriptions. For further information, please read chapter 9.3 .
2. Then press the button on the right-hand side Save template


Figure 9.3.4.: Save Template pop-up menu

3. Enter a unique template name and confirm this with OK .

9.3.2. Load Template

The device has the option to create and load templates/templates for corresponding measurement sequences. To create a template, please read chapter 9.3.1. To load an existing template, please follow the steps below:

1. In sample mode, press the button Load template .



Figure 9.3.5.: Load template menu

2. Select the desired template from the list displayed and confirm with OK .
3. The template is loaded and is now ready for measurement.

9.4. Qualification

This menu contains all qualification options for the device. The following options are available:


Figure 9.4.1.: Menu option selection Qualification

Options	Description
TOC	TOC calibration with sucrose
SST	System suitability test with sucrose and benzoquinone
Conductivity	Calibration of the conductivity sensors

9.4.1. Calibration TOC

To achieve a sufficient operating temperature, the system must be switched on one hour before starting the calibration. The calibration solutions must be warmed up to at least room temperature beforehand. If temperatures below 18°C are displayed at the start of the measurements, the measurement must be interrupted until the solutions have reached at least room temperature. The optimum temperature range is 20°C to 25°C as the input temperature for calibration.

To calibrate the TOC value, all current measurements must be stopped. In the Qualification selection menu, press **TOC**, see illustration 9.4.1


Figure 9.4.2.: User interface for TOC calibration

Number	Display	Description
1	Table field	Displays all measurements, including concentration [ppb], type, port, conductivity [µS/cm] and time.
2	Result field	Displays the result of the measurement, including factor, offset, coefficient of determination R2 and the percentage deviation from the current value.
3	Info field	Displays information on the current measurement.

The table shows the required solutions in the corresponding concentrations. Each solution is measured five times as standard. The measurement duration and measurement repetitions can be adjusted in the settings, please refer to chapter 9.4.2.

You can use the info field to view useful information about the current measurement.

- Current Conc.(ppb) - Solution concentration of the current measurement sample
- Next Step - Next measurement step
- Status - Device status [Idle = idle, Flush = flush, measurement = measurement]
- Elapsed-Time - Displays the elapsed time
- User - Displays the logged in user
- ID - displays the measurement ID of the current measurement

To start the measurement, check that the solutions are connected to the correct port and that the connection is hand-tight. A loose connection can falsify the measurement results.

Start the calibration of the device by pressing the **Start** -buttons.

Once the measurement is complete, the result is displayed on the screen, Figure 9.4.3.


Figure 9.4.3.: Result of a successfully completed calibration

Number	Designation	Description
1	Factor	Factor for TOC calibration
2	Offset	Calibration offset
3	R ² (coefficient of determination)	Goodness of fit

Via the Update factor? the new calibration factor is confirmed. The pop-up window shows the new calibration factor, the previous calibration factor and the percentage deviation between the values.

With the -button to cancel the calibration process. The previous calibration factor is retained.

9.4.2. Calibration-Sequence

The device offers the option of adjusting the calibration as required. To use this function, navigate to the main menu and click on . Then click on the button . The following user interface opens:

Calibration Sequence

Conc.(ppb)	Type	Input
1000	Sucrose	Input 3
0	Water	Input 4
NA	Water	Input 1
NA	Water	Input 1
NA	Water	Input 1
NA	Water	Input 1

Iterations

5

OK

Cancel

Figure 9.4.4.: Clibration sequence menu

The concentration, type and port can be changed via the individual drop-down lists. The following settings are provided as standard for the calibration:

→ 1000 ppb - Sucrose - Port 1

→ 0 ppb - Water - Port 2

If you want to change the concentration, please click on the corresponding drop-down menu "Conc. (ppb)".

Please note, however, that this shifts the upper detection limit and the device is then qualified for the set concentration range.

To add further calibration points, click on "NA". You can use this function to perform a multi-point calibration. Please note that you must provide the set solutions with the corresponding concentrations for calibration.

Calibration Sequence

Conc.(ppb)	Type	Input
1000	Sucrose	Input 3
0	Water	Input 4
100	Water	Input 1
200	Water	Input 1
300	Water	Input 1
400	Water	Input 1
500		
600		
700		
800		
900		
✓ 1000		
1100		
1200		
1300		

5

Cancel

Figure 9.4.5.: Drop-down selection list Calibration sequence

If you want to change the type of substance for the calibration, click on the "Type" drop-down menu.

Please note, however, that some substances are more difficult to break down and therefore change the measured values.

Calibration Sequence

Conc.(ppb)	Type	Input
1000	Sucrose	Input 3
0	Water	Input 4
NA	Water	Input 1
NA	Water	Input 1
NA	Water	Input 1
NA	Water	Input 1

Iterations

5

OK

Cancel

Figure 9.4.6.: Drop-down selection list Calibration sequence

The port of the solution can be set via the "Port" drop-down menu. During calibration, make sure that all solutions are connected to the correct ports. The device automatically pulls the calibration solution into the system.

Calibration Sequence

Conc.(ppb)	Type	Input
1000	Sucrose	Input 3
0	Water	Input 1
NA	Water	Input 2
NA	Water	✓ Input 3
NA	Water	Input 4
NA	Water	Input 1

Iterations

5

OK

Cancel

Figure 9.4.7.: Drop-down selection list Calibration sequence

If the number of calibration points is greater than the number of ports in the device, the device automatically waits and informs the user that a solution change is required.

9.4.3. SST

To achieve a sufficient operating temperature, the system must be switched on one hour before the start of the SST. The test solutions must be warmed up to at least room temperature beforehand. If temperatures below 18°C are displayed at the start of the measurements, the measurement must be interrupted until the solutions have reached at least room temperature. The optimum temperature range is 20°C to 25°C as the input temperature for the SST.

Please only use solutions for the SST in accordance with your local and national standards. Information on the individual solutions can be found in your national pharmacopoeias and guidelines.

The device has the option of carrying out a system suitability test (SST for short). To do this, all current measurements must be interrupted. In the selection menu, press the button **SST** (see illustration 9.4.1 in chapter 9.4


Figure 9.4.8.: SST user interface

Number	Display	Description
1	Table field	Displays the values of all measurements, with concentration [ppb], type, port, measured concentration [ppb] and time.
2	Result field	Displays the result of the measurement, including TOC values for ultrapure water, sucrose, 1,4-benzoquinone and the percentage result of the SST.
3	Info field	Displays information on the current measurement.

A pop-up window prompts you to provide the solutions, Figure 9.4.9. Connect them to the device as specified and tighten the connections hand-tight. Confirm the provision of the solutions with .


Figure 9.4.9.: Pop-up window prompts you to provide the solutions

The table shows the measurements with port and expected concentration of the solution. Five measurements are carried out as standard. The number of measurements can be adjusted in the settings, please refer to chapter ***.

Start the measurement by pressing the -buttons. Once the automatic measurement process is complete, a pop-up window appears indicating that the measurement has been completed. After confirmation, the results are displayed in the graphical user interface, Figure 9.4.10.


Figure 9.4.10.: Results of the SST (System Suitability Test)

9.4.4. Calibration of conductivity sensors cell constant

To achieve a sufficient operating temperature, the system must be switched on one hour before starting the calibration. A TOC measurement or operation of the device with a UV lamp must be interrupted one hour before the start of calibration in order to avoid temperature changes in the sensors during conductivity calibration. The calibration solutions must be warmed up to at least room temperature beforehand. If temperatures below 18°C are displayed at the start of the measurements, the measurement must be interrupted until the solutions have reached at least room temperature. The optimum temperature range is 20°C to 25°C as the input temperature for calibration.

To calibrate the cell constant of the conductivity sensors, all current measurements must be stopped. In the Qualification selection menu, press **Conductivity**, see illustration 9.4.1


Figure 9.4.11.: User interface for calibrating the cell constant of the conductivity sensors

Number	Display	Description
1	Table field	Displays all measurements, including conductivity of standard [$\mu\text{S}/\text{cm}$], type, port, measured conductivity [$\mu\text{S}/\text{cm}$] and time.
2	Result field	Displays the result of the measurement, including factor and the percentage deviation from the current value.
3	Info field	Displays information on the current measurement.

The table shows the required and measured conductivity of the standard. Each solution is measured five times as standard. The measurement duration and measurement repetitions can be adjusted in the settings, please refer to chapter 10.1.

You can use the info field to view useful information about the current measurement.

- Current Cond.($\mu\text{S}/\text{cm}$) - Conductivity of the standard as measurement sample
- Status - Device status [Idle = idle, Flush = flush, measurement = measurement]
- Elapsed-Time - Displays the elapsed time
- User - Displays the logged in user
- ID - displays the measurement ID of the current measurement

To start the measurement, check that the solutions are connected to the correct port and that the connection is hand-tight. A loose connection can falsify the measurement results.

Start the calibration of the device by pressing the **Start**-buttons.

The message window indicates the conductivity of the standard to be connected, Figure 9.4.12 . The conductivity of the connected standard at 25 °C must match the value displayed in the message window.


Figure 9.4.12.: Message window at the beginning of the calibration of the factors for the cell constants

Once the measurement is complete, the result is displayed on the screen, Figure 9.4.13


Figure 9.4.13.: Result of a successfully completed calibration of the cell constant factors

Nummer	Bezeichnung	Beschreibung
1	Old-C1-Factor	previous factor for the calibration of the factors for the cell constant of conductivity sensor 1
2	New-C1-Factor	new factor for the calibration of the factors for the cell constant of conductivity sensor 1
3	Correction	Deviation of factor C1 from the previous value
4	Old-C2-Factor	previous factor for the calibration of the factors for the cell constant of conductivity sensor 2
5	New-C2-Factor	new factor for the calibration of the factors for the cell constant of conductivity sensor 2
6	Correction	Deviation of factor C2 from the previous value

Via the Update factor? the new calibration factor is confirmed. The pop-up window shows the new calibration factor, the previous calibration factor and the percentage deviation between the values.

With the -button to cancel the calibration process. The previous calibration factor is retained.

9.5. Direct

Direct mode is the manual mode of the device. In this measuring mode, the following components can be operated/adjusted by pressing buttons:

- Port 1 - 4
- Valves 1 - 3
- UV lamp
- Pump 1 - 2 (ON/OFF, H: high/low pump speed)
- IO outputs 1 - 4

To access Direct mode, navigate to the main menu and click on the button **Direct**. The following user interface then opens:


Figure 9.5.1.: Direct user interface menu

Number	Description
1	Switching port 1 - 4
2	Switching valves 1 - 3
3	ON/OFF UV lamp
4	Sensor status Leakage in the device
5	ON/OFF Pump 1 - 2, H: high, low pump speed
6	Conductivity value, temperature compensation method and temperature value in conductivity sensor 1 and 2
7	Conductivity difference C2-C1
8	Instantaneous values for pressure and system temperature, with selected option
9	Default values for IO outputs 1 - 4
10	Operating options for the autosampler, with selected option
11	Switches to the graphical display of the measurement process
12	Switches to the Quicksetup user interface, chapter 10.1
13	Returns to the start page

With a click on the button **Graph** to access the chart view, shown in Figure 9.5.2 . There you can see the course of the conductivity values of the two sensors C1 and C2, as well as the temperature of the two sensors.


Figure 9.5.2.: Diagram of the measurement process

The time displayed on the X-axis can be adjusted using the selection button below the diagram. The following options are available: 10 min, 1h (1 hour), 6h (6 hours), 12h (12 hours). Via the button Direct you get back to the "Direct" user interface.

9.6. About

Use the button **About** to display general device information. The button **About** is available on the main menu.


Figure 9.6.1.: Menu About

Info	Description
SN	Serial number of the device
Version	Software version of the device
MAC address	MAC address of the device
OS	Operating system
Firmware	Controller firmware
Computer name	Computer-Name
IP Address	IP-Adress of the device
HW config	Hardware configuration code
HW Rev	Hardware revision-number

9.7. Status

General information on the status of the device can be displayed in the status overview. This includes operating hours of individual sensors, pumps and the UV lamp. To access the status overview, navigate to the main menu. Then click on the button **Status**. The graphical user interface is then displayed.


Figure 9.7.1.: Status menu


Number	Description
1	Running time of the UV lamp with display as a pie chart and date of the last change
2	Display of device errors as error codes
3	Device information (serial number, version, etc.)
4	Running times of the pump hoses of pumps 1 and 2 with display as a pie chart and date of last replacement
5	Running times of all valves, pumps and UV-EVG with display as progress bar and date of last replacement
6	Button OK Returns to the main menu

9.8. Measurement data backup and visualization

The device offers various options for exporting and analyzing data. You can export data cloud-based or to an external removable storage device (chapter: 9.8.2) or have the measurement data output as a PDF report (chapter: 9.8.3).

9.8.1. Cloud-based export of data

Cloud export allows the user to provide data for the service and automatically transfer it to the Membrapure Cloud. To do this, navigate to the main menu and click on [More](#). Then click on [Cloud](#). The user interface opens.

The screenshot displays the 'Cloud menu' interface. At the top, there is a 'Select Preset' dropdown menu currently set to 'All'. Below this, there are two date selection fields: 'From' (21.05.2025 00:00:00) and 'To' (23.05.2025 23:59:59). A duration field shows '03 day(s)'. Three buttons are available: 'Today', 'Last 3d', and 'Last 7d'. Under the 'Path local' section, a text field contains the path 'C:\TOC...\a08t0001-all-20250523-120004'. The 'Online path name' section has a text field with 'a08t0001-all-20250523-120004'. At the bottom, there are two large blue buttons: 'Export' (with a cloud icon) and 'Home' (with a house icon).

Figure 9.8.1.: Cloud menu

The following templates (presets) are available:

Preset	Description
Qualification	Exports only data created during the qualification of a device
Measurement	Exports only data that was recorded as a measurement
All	Exports the qualification and measurement data
Service	In addition to the measurement and qualification data, debug data is also transmitted if it has been recorded using the corresponding settings

The "From" and "To" buttons can be used to define the period to be exported or transferred.

Button	Description
From	A date selection is opened that defines the day, month and year. This represents the start date of the export
To	A date selection opens, which defines the day, month and year. This represents the end date of the export

The small display below shows the set number of days.

Below the settings, three presets can be used to preset a period that starts with the current date and ends with the period indicated on the buttons.

The Path display indicates the local path on the device under which the data is buffered before being exported to the cloud.

Once all settings have been set correctly, you can use the **Export** to start the transfer. Once this has been successfully completed, a brief confirmation appears. The button **Home** cancels the process and the display returns to the main display.

9.8.2. Export measurement data

The export allows the user to export data for a specific period and make it available to the service if required. There are various options for this under the "Presets" selection, which determine the scope of the transmitted data. The data is written to an external data carrier. The target directory is the "Export" folder by default. Connect an external removable storage device to the device. The USB 2.0 port on the right-hand side of the device is used for this purpose (see chapter) ?? The following templates (presets) are available:

Preset	Description
Qualification	Exports only data created during the qualification of a device
Measurement	Exports only data that was recorded as a measurement
Sample	Exports data that was generated as a single measurement under Sample
All	Exports the qualification and measurement data

Service	In addition to the measurement and qualification data, debug data is also transmitted if it has been recorded using the corresponding settings
---------	--

The "From" and "To" buttons can be used to define the period to be exported or transferred.

Button	Description
From	A date selection is opened that defines the day, month and year. This represents the start date of the export
To	A date selection opens, which defines the day, month and year. This represents the end date of the export

The small display below shows the set number of days. Below the settings, three presets can be used to preset a period that starts with the current date and ends with the period indicated on the buttons. With "Select Path" an external data carrier is selected to which the data is written after export. The Path display indicates the local path on the data carrier under which the data will be written after export.

Once all settings have been set correctly, you can use the **Export** to start the transfer. Once this has been successfully completed, a brief confirmation appears. The button **Home** cancels the process and the display returns to the main display.

9.8.3. Display and report output of data

The data can also be viewed in the Data Viewer, chapter 9.9 and output as a preset log in PDF format. To use this function, please navigate to the main menu, see chapter 7. Click here on the button **Data Viewer**.

9.9. DataViewer

The DataViewer can be used for several purposes. Reports can be created, measurements checked or error reports analyzed. Depending on the objective, you can navigate to the DataViewer modules. Start with the data input dialog **Data** Chapter 9.9.1.


Figure 9.9.1.: Menu buttons of the DataViewer

Number	Button	Description
1	Data	Shows the data entry dialog
2	Live	Shows the data display dialog
3	Preview	Shows the document preview
4	Export	Shows the data export dialog
5	Home	Exits the Data Viewer module

9.9.1. Data input dialogue

In the data entry dialog **Data** the templates, time period and other options are first selected before further steps can be taken, e.g. to create a report. If no data is available for report creation in the selected time period, the upper menu buttons remain grayed out.


Figure 9.9.2.: Data entry dialog

The first step is to select a template (preset). Each preset refers to a data type, defines which graphs are to be displayed, which data is to be calculated and specifies the design of the reports. A distinction is made between the data types:

- Qualification measurement Calibration
- Qualification measurement System Suitability Test (SST)
- Monitoring
- Log messages Log
- Error log messages ErrorLog

Template name	Description	Measurement run selection
Calibration	Results of a calibration measurement. The results can be seen under "Live" -> "Info" and in the report they are on one page	Individual selection
Calibration Accepted	Like the "Calibration" template, except that only successful calibrations are listed under the measurement run selection	Individual selection
SST	Results of an SST measurement. The results can be seen under "Live" -> "Info" and in the report they are on one page	Individual selection
SST Passed	Like the "SST" template, except that only successful SSTs are listed under the measurement run selection	Individual selection
Monitor	The data from continuous measurements is displayed here. All (standard) or several measurements can be selected in the measurement run selection.	Multiple selection
Log	All log messages of the device in the selected period are displayed here	invisible
Log Error	All error messages of the device in the selected period are displayed here	invisible
Log User	All messages from the device during the selected period that are generated by explicit operation by the user are displayed here	invisible
Log Archive	As with "Log", all log messages from the device are displayed here. However, the lines in the report are as narrow as possible with a small font so that as many lines as possible fit on one page.	invisible

Only for qualification or monitor measurements can and must the measurement runs to be examined be selected under "Selected Measurements", otherwise the area is invisible.

In the data entry dialog, the selected period for the data selection is displayed and defined under "From - To", chapter 9.9.2. This period also determines which measurement runs are displayed in the measurement run selection below.

If the measurement run selection ("Selected Measurements") is visible, all measurement runs for the selected time period are displayed. Depending on the template selection, at least one run can and must be selected. In the case of templates such as Monitor, several or all can also be selected. You can change the selection by clicking on the respective run. The status fields at the bottom left show whether and how much data could be extracted from the database with the selected template, date range and measurement runs, how much progress has been made with their processing and, depending on the template, possibly also other display data. If the field is green, this means that the selection generates data to the permitted extent or that the data has been fully processed. If the field is red, too little or too much data has been selected or the calculation has not yet been completed.

Field	Description	Marked in red if
Row Count	Number of result rows. If no data could be extracted - i.e. the field shows 0 rows - the line is red, then all other menu fields at the top are grayed out. If too much data has been extracted, the field is also red. The menu buttons at the top are NOT greyed out and the areas of the program are available. However, it should be noted that the data is incomplete in this case. The limit for the maximum permitted number of data can be set in the main program.	No or too much data
Page Count	Number of result pages in the report. If no data could be extracted - i.e. the field shows 0 pages - the line is red, then all other menu fields at the top are grayed out. If too much data was extracted, the field is also red. The menu buttons at the top are NOT greyed out and the areas of the program are available. However, it should be noted that the data is incomplete in this case. The limit for the maximum permitted number of data can be set in the main program.	no or too many pages
Preview	Displays the calculation progress of the report preview. As long as this is not set to 100%, the above menu fields remain grayed out.	less than 100%
Live	Displays the calculation progress of the graphs that are visible in both the data display (Live) and the report. As long as this is not set to 100%, the above menu fields remain grayed out.	less than 100%
Export	Displays the calculation progress of the complete report. As long as this is not at 100%, you can access the export dialog, but not yet export.	less than 100%

9.9.2. Date selection dialogue

If you press the "From - To" period, the date selection dialog opens, Figure 9.9.3in which a new date range can be set.


Figure 9.9.3.: Date selection dialog

Number	Description
1	Start date of the period
2	End date of the period
3	Buttons for changing the date, depending on whether the start or end date is selected
4	"OK" confirms the selection; "Cancel" exits the date selection menu without making any changes

First select what you want to change - the start date or the end date. By pressing buttons, area [3] in figure 9.9.3 changes the selected date or both:

Designation	Description	Concerns
1 day	one day more	Selection
-1 day	one day less	Selection
1 month	one month more	Selection
-1 month	one month less	Selection
1 year	one year more	Selection
-1 year	one year less	Selection
Last 3 days	the last 3 days until today	Both
Last week	the last week until today	Both
Last month	the last month until today	Both

9.9.3. Data display dialogue

You don't always want to create a report from the data, but rather view it directly. Use the data display dialog **Live** for this purpose. All data from the selected area is displayed here with the selected default settings - in table form or as a graph. Information calculated from the data is also listed in the "Info" area.


Figure 9.9.4.: Data display dialog in the table view

If there are no graphs in a template, e.g. in the log messages (Log), the Chart button is grayed out. To the right of the "Info" button, you can select the table or graph to be displayed, if there are several, Figure 9.9.5.


Figure 9.9.5.: Selection for multiple charts and tables

The views are available:

View	Description
Table view	In the table view (Table), the tables that would also appear in the report can be viewed directly. The columns are therefore not static, but vary depending on the table. Not all rows are always loaded at the beginning, but only when you scroll to the last row.
Chat view	The "Chart" graph view shows all graphs that would also appear in the report.
Calculated or individual information	Under the "Info" view, all calculated or individual information determined by the template is visible. This can be measurement results, minima or maxima or even the size of the database, Figure 9.9.7 . Depending on the template, some of this information can also be found in the report.


Figure 9.9.6.: Data display dialog in the diagram view

The graph is interactive and can be moved in all four directions by swiping. To the right are five additional control buttons that can be used to zoom in and out of the graph or return to the original view.

Field	Description
	Enlarging the view along the X-axis
	Reducing the view along the X-axis
	Reset the view to the entire graph
	Enlarging the view along the Y-axis
	Reducing the view along the Y-axis


Figure 9.9.7.: Data display of calculated or individual information

9.9.4. Document preview

In the document preview, tab **Preview**, a preview of the document is displayed as it was previously configured in the data entry dialog.

If the document contains chart objects, these can be displayed beforehand via the data display dialog in the chart view, tab **Chart**. The report to be created can be viewed in the report preview. In order to be able to create the preview as quickly as possible, only a maximum of two pages per chapter can be viewed here by default.


Figure 9.9.8.: Document preview

9.9.5. Data export dialogue


Figure 9.9.9.: Data export dialog

At this point, once the complete report has been calculated, it can be exported or printed in various formats. If the report has not yet been calculated, all buttons are grayed out and show the calculation progress. Even if a report is exported, the corresponding button is grayed out until the process is complete. Data can be exported as HTML, CSV or PDF. HTML and PDF are formats in which the report is output in a paper layout that looks exactly like a printout. HTML outputs a folder with all content and is the much faster of the two exports. CSV outputs a comma-separated text file, which in turn can be imported into programs such as MS Excel to create your own calculations or graphs. When printing (Print), the default printer preset in the service area is used with the default settings set there.


Figure 9.9.10.: Drive selection dialog

After pressing one of the three export buttons in the report export window, the small drive selection window appears. Here you can select the drive to which you want to export. The result will then be found in this drive in the "Export" folder. The name of the result file or folder contains both the creation date and the template name.

9.10. Remote access

The device is able to establish a remote session via an existing Internet connection. To use this function, navigate to the main menu, see chapter 7 .

There please click on the button Remote . An AnyDesk remote session opens.

9.11. Messages

The device can display various error messages in the main menu.


Figure 9.11.1.: Info menu Messages

Each error is displayed with an error code and a brief description. In service cases, please make a note of the error code and pass it on to the relevant service employee.

Error code	Meaning
------------	---------

Once the error has been corrected, you can view the message by clicking on the **Message** button. Confirm the pop-up window with **OK**.

9.12. Shut down system

To shut down the system, please navigate to the main menu. There you will find the button **Shutdown**.


Figure 9.12.1.: Shutdown button

By pressing the button **Shutdown** ensures a safe shutdown of the system. If the system is shut down, turn the On-Off cam switch to the "0" position.

Chapter 10

Settings

10.1. Quick setup

The user can change various device settings in this menu. To access the Quick setup menu, please navigate to the main menu and click on the button Quick setup . Then log in with your user name and password. The graphical user interface of the quick setup menu opens. Various option tabs can be selected via the tabs.

License	Sensor Signal Processing I	Sensor Signal Processing II	IO Details
User Interface	System	Operation Limits (Warning)	Operation Limits (Alarm)

Figure 10.1.1.: Quicksetting menu

The following settings can be made via the **User interface** set:

Setting	Description
Default Language	Setting the device language
Accuracy TOC Value	Number of decimal places for the display of the TOC value
Name of Input 1	Name for the fluid input at port 1. This name is then used for the various measurement modes.
Name of Input 2	Name for the fluid input at port 2. This name is then used for the various measurement modes.
Name of Input 3	Name for the fluid input at port 3. This name is then used for the various measurement modes.
Name of Input 4	Name for the fluid input at port 4. This name is then used for the various measurement modes.
Activate sound?	Switches acoustic notification for certain events on or off
Enable autostart?	If the autostart option is switched on, the system automatically returns to monitor mode after a restart and continues the measurement
Show conductivity in μS ?	If this option is activated, the measured values for conductivity are displayed in the unit $\mu\text{S}/\text{cm}$


Figure 10.1.2.: "User Interface" tab in the Quicksetup menu

Via the tab **System** it is possible to make system settings on the device.

Setting	Description
Enable touchscreen keyboard?	Switches the virtual keyboard on or off
Auto Logout time	Time in s until the logged-in user is automatically logged out
Monitor mode recording interval	Time interval between two measurements in monitor mode

The limit values for warning messages for the following parameters are set on the **Operation Limits (Warning)** setup:

Setting	Description
TOC1 warning level	Limit value from which a warning is issued for TOC1; limit value for the warning must not exceed the alarm limit value
TOC2 warning level	Limit value from which a warning is issued for TOC2; limit value for the warning must not exceed the alarm limit value
TOC3 warning level	Limit value from which a warning is issued for TOC3; limit value for the warning must not exceed the alarm limit value
Temperature1 warning level	Limit value above which a warning is issued for the temperature in the conductivity sensor1 ; the limit value for the warning must not exceed the alarm limit value

Temperature2 warning level	Limit value above which a warning is issued for the temperature in the conductivity sensor2 ; the limit value for the warning must not exceed the alarm limit value
Conductivity1 warning level	Limit value above which a warning is issued for the conductivity measured value in the conductivity sensor1 ; the limit value for the warning must not exceed the alarm limit value
Conductivity2 warning level	Limit value above which a warning is issued for the measured conductivity value in the conductivity sensor2 ; the limit value for the warning must not exceed the alarm limit value
Pressure sensor warning level	Limit value above which a warning is issued for the inlet pressure measurement; the limit value for the warning must not exceed the alarm limit value

There are certain value ranges in which the limit values for warning messages can be set. These are specified in the settings. The limit values for alarm messages for the following parameters are set on the **Operation Limits Alarm** setup:

Setting	Description
TOC1 alarm level	Limit value above which an alarm is issued for TOC1; the alarm limit value must not be below the limit value for the warning.
TOC2 alarm level	Limit value above which an alarm is issued for TOC2; the alarm limit value must not be below the limit value for the warning.
TOC3 alarm level	Limit value above which an alarm is issued for TOC3; the alarm limit value must not be below the limit value for the warning.
Temperature1 alarm level	Limit value above which an alarm is issued for the temperature in the conductivity sensor1 and the measurement is stopped; the alarm limit value must not be below the limit value for the warning.
Temperature2 alarm level	Limit value above which an alarm is issued for the temperature in the conductivity sensor2 and the measurement is stopped; the alarm limit value must not be below the limit value for the warning.
Conductivity1 alarm level	Limit value above which an alarm is issued for the conductivity in the conductivity sensor1; the alarm limit value must not be below the limit value for the warning.
Conductivity2 alarm level	Limit value above which an alarm is issued for the conductivity in the conductivity sensor2; the alarm limit value must not be below the limit value for the warning.
Pressure sensor alarm level	Limit value above which an alarm is issued for the inlet pressure measurement; the limit value for the warning must not exceed the alarm limit value

There are certain value ranges in which the limit values for the alarms can be set. These are specified in the settings.

10.2. Device settings

10.2.1. Change UV

This function should only be carried out if the UVC lamp has been replaced. If this function is carried out without replacing the lamp, the device may no longer be able to check the current lamp correctly.

Instructions for removing and replacing the UV lamp can be found in chapter 13.1. This function is used to reset the operating hours of the UV lamp. If you want to reset the operating hours, confirm the prompt with .


Figure 10.2.1.: Resetting the operating hours of the UV lamp

10.2.2. Change Pump Tube

This function should only be carried out if the pump hose has been replaced. If this function is carried out without replacing the pump hose, the device may no longer be able to check the current pump hose correctly.

Instructions for removing and replacing the pump hose can be found in chapter 13.2.1. This function is used to reset the operating hours of the pump hose. If you want to reset the operating hours, confirm the prompt with .


Figure 10.2.2.: Resetting the operating hours of pump 1


Figure 10.2.3.: Resetting the operating hours of pump 2

10.2.3. Backup

The device is able to create backups to ensure greater data security.

To create a backup, navigate to the main menu and click on the button More. Then click on Backup.


Figure 10.2.4.: Backup pop-up window

The device creates a backup and informs the user via a pop-up window with the file name and file path.

10.2.4. Update Software

If a new software version is installed on the device, the data and settings remain unaffected.

It is recommended to carry out an automatic backup before each software update. To carry out a backup, please read chapter 10.2.4 and follow the instructions described therein.

To carry out a software update, navigate to the main menu. There, click on the button **More**. Then press the button **Update**. Confirm the pop-up window with **Yes**.


Figure 10.2.5.: Software update pop-up window

Then select the update file from your removable storage device (USB stick) and press **OK**.


Figure 10.2.6.: Selecting the update file

The device updates automatically and starts the new software automatically. Check that the installation was successful by clicking on **About** and check the serial number.

10.2.5. License

There are several ways to purchase software extensions. For an up-to-date list, please visit the website <https://membrapure.de/de/> or see the chapter 14 .

To enter a purchased license key and activate the associated software, please navigate to the main menu. Click on the button More and then click on the button License .


Figure 10.2.7.: License menu

Enter the license key and confirm the entry with OK .

Chapter 11

Factory settings

The device offers the option of resetting various parts of the device to the factory settings. These include

- Reset user
- Reset rights
- Reset settings
- Reset device (includes all methods mentioned above)

These functions can cause considerable damage to the device. All deleted data is irrevocably lost. For this purpose, you must go through a 2-factor authentication process to confirm.

11.1. User reset

To delete all set users and restore the factory settings, use the "User Reset" function. To use this function, press the button in the main menu **More**. Then press the button **Factory Reset User**.


Figure 11.1.1.: Pop-up window User factory settings

If you really want to delete the user data, confirm the request with **Yes**. A new pop-up menu then appears with a "Code" entry. Make a note of the number displayed and contact the company's service department membraPure GmbH. Further information on 2-factor authentication can be found in the chapter 12.2. After you have confirmed your code entry, the settings are reset. To create a new user, please read chapter 8.3.

11.2. Rights reset

To delete all set rights and restore the factory settings, use the "Rights Reset" function. To use this function, press the button in the main menu **More**. Then press the button **Factory reset rights**.


Figure 11.2.1.: Pop-up window Rights factory settings

If you really want to delete the current rights management, confirm the prompt with **Yes**. A new pop-up menu then appears with a "Code" entry. Make a note of the number displayed and contact the company's service department membraPure GmbH. Further information on 2-factor authentication can be found in the chapter 12.2. After you have confirmed your code entry, the settings are reset. To set new rights, please read chapter 8.2.

11.3. Settings reset

To delete all existing settings and restore the factory settings, use the "Factory reset settings" function.

To use this function, press the button in the main menu **More**. Then press the button **Factory reset settings**.


Figure 11.3.1.: Pop-up window Device settings Factory settings

If you really want to delete the current settings, confirm the prompt with **Yes**. A new pop-up menu then appears with a "Code" entry. Make a note of the number displayed and contact the company's service department membraPure GmbH . Further information on 2-factor authentication can be found in the chapter 12.2. After you have confirmed your code entry, the settings are reset. To make new settings on the device, please read chapter 10.1 and chapter 10.2.

11.4. Device factory settings

This function includes resetting the following device areas:

- Rights management
- User administration
- Device settings

This function resets the device to the factory settings. To use this function, press the button in the main menu **More** . Then press the button **Factory reset all** .


Figure 11.4.1.: Pop-up window Devices factory settings

If you want to reset the device to the factory settings, confirm the prompt with . A new pop-up menu then appears with a "Code" entry. Make a note of the number displayed and contact the company's service department membraPure GmbH . Further information on 2-factor authentication can be found in the chapter 12.2.

After you have confirmed your code entry, the device is reset to the factory settings. To set up the device again, please read chapter 8 .

Chapter 12

Repair and service

If you have any questions or problems that are not described in this manual, please contact our service department.

Address:

membraPure GmbH
Wolfgang-Küntschers-Straße 16
16761 Hennigsdorf
Germany

Service telephone: +49 3302 - 201 20 0

12.1. Cloud

To upload measurement data to the cloud, navigate to the main menu and click on **More**. Then click on **Cloud**. You can use the user interface as described in chapter 9.8.1 for cloud based transfer of data. Our service department can use this data to initiate the necessary troubleshooting measures.

12.2. 2-factor authentication

2-factor authentication is required to use some of the device's functions. This protects the user against unwanted changes to the device. This includes, for example

- User reset (see chapter 11.1)
- Rights reset (see chapter 11.2)
- Settings reset (see chapter 11.3)
- Device factory setting (see chapter 11.4)


Figure 12.2.1.: Pop-up window 2-factor authentication

You will be asked to enter a code for the secure functions of the device. First contact the service team of the company membraPure GmbH .

Then tell the employee the code displayed on the device. You will then receive a new code, which you enter in the bottom field of the pop-up menu.

Confirm the entry with .

12.3. Error numbers and displays

Any errors that occur are registered by the software and assigned error numbers. The error numbers appear on the main screen, in the status information and are exported with the data. The table shows the list of occurring error numbers with a description:

Error code	Category	Level	Show-dialog	Show status	ignore
err_10001		yellow	false	true	false
err_10002		yellow	false	true	false
err_10003		yellow	false	true	false
err_10004		yellow	false	true	false
err_10005		yellow	false	true	false
err_10006	system_error	yellow	false	true	false
err_10007	system_error	yellow	false	true	false
err_10008	system_error	yellow	false	true	false
err_10009	system_error	yellow	false	true	false
err_10010	system_error	yellow	false	true	false
err_10011	system_error	yellow	false	true	false
err_10012	system_error	yellow	false	true	false
err_10013	system_error	yellow	false	true	false
err_10014	system_error	yellow	false	true	false
err_10015	system_error	yellow	false	true	false
err_10016	system_error	yellow	false	true	false
err_10017	system_error	yellow	false	true	false
err_10018	system_error	yellow	false	true	false
err_10019	system_error	yellow	false	true	false
err_10020	system_error	yellow	false	true	false
err_10021	system_error	yellow	false	true	false
err_10022	system_error	yellow	false	true	false
err_10023	system_error	yellow	false	true	false
err_10024	system_error	yellow	false	true	false
err_10025	system_error	yellow	false	true	false
err_10026	process_error	yellow	false	true	false
err_10027	system_error	yellow	false	true	false
err_10028	system_error	yellow	false	true	false
err_10029	system_error	yellow	false	true	false
err_10030	system_error	yellow	false	true	false

Error code	Category	Level	Show-dialog	Show status	ignore
err_10031	system_error	yellow	false	true	false
err_10032	system_error	yellow	false	true	false
err_10033	system_error	yellow	false	true	false
err_10034	system_error	yellow	false	true	false
err_10035	system_error	yellow	false	true	false
err_10036	system_error	yellow	false	true	false
err_10037	system_error	yellow	false	true	false
err_10038	system_error	yellow	false	true	false
err_10039	process_error	yellow	false	true	false
err_10040	process_error	yellow	false	true	false
err_10041	process_error	yellow	false	true	false
err_10042	process_error	yellow	false	true	false
err_10043	process_error	yellow	false	true	false
err_10044	process_error	yellow	false	true	false
err_10045	process_error	yellow	false	true	false
err_10046	process_error	yellow	false	true	false
err_10047	process_error	yellow	false	true	false
err_10048	process_error	yellow	false	true	false
err_10049	process_error	yellow	false	true	false
err_10050	process_error	yellow	false	true	false
err_10051	process_error	yellow	false	true	false
err_10052	process_error	yellow	false	true	false
err_10053	process_error	yellow	false	true	false
err_10054	process_error	yellow	false	true	false
err_10055	process_error	yellow	false	true	false
err_10056	process_error	yellow	false	true	false
err_10057	process_error	yellow	false	true	false
err_10058	process_error	yellow	false	true	false
err_10059	process_error	yellow	false	true	false
err_10060	process_error	yellow	false	true	false
err_10061	process_error	yellow	false	true	false
err_10062	process_error	yellow	false	true	false
err_10063	process_error	yellow	false	true	false
err_10064	system_error	yellow	false	true	false
err_10065	process_error	yellow	false	true	false
err_20001	system_error	red	true	true	false
err_20002	system_error	red	true	true	false
err_20003	system_error	red	true	true	false
err_20004	process_error	red	true	true	false
err_20005	system_error	red	true	true	false

Error code	Category	Level	Show-dialog	Show status	ignore
err_20006	process_error	red	true	true	false
err_20007	process_error	red	true	true	false
err_20008	system_error	red	true	true	false
err_20009	process_error	red	true	true	false
err_20010	process_error	red	true	true	false
err_20011	process_error	red	true	true	false
err_19001	board_error	yellow	false	true	false
err_19002	board_error	yellow	false	true	false
err_19003	board_error	yellow	false	true	false
err_19004	board_error	yellow	false	true	false
err_19005	board_error	yellow	false	true	false
err_29006	board_error	red	true	true	false
err_29007	board_error	red	true	true	false
err_29008	board_error	red	true	true	false
err_29009	board_error	red	true	true	false
err_19010	board_error	yellow	false	true	false
err_29011	board_error	red	true	true	false
err_29012	board_error	red	true	true	false
err_29013	board_error	red	true	true	false
err_29014	board_error	red	true	true	false
err_29015	board_error	red	true	true	false
err_19016	board_error	yellow	false	true	false
err_19017	board_error	yellow	false	true	false
err_19018	board_error	yellow	false	true	false
err_19019	board_error	yellow	false	true	false
err_29020	board_error	red	true	true	false
err_29021	board_error	red	true	true	false
err_29022	board_error	red	true	true	false
err_29023	board_error	red	true	true	false
err_29024	board_error	red	true	true	false
err_29025	board_error	red	true	true	false
err_19026	board_error	yellow	false	true	false
err_19027	board_error	yellow	false	true	false
err_19028	board_error	yellow	false	true	false
err_19029	board_error	yellow	false	true	false
err_29030	board_error	red	true	true	false
err_29031	board_error	red	true	true	false
err_29032	board_error	red	true	true	true
err_29033	board_error	red	true	true	false
err_29034	board_error	red	true	true	false

Error code	Category	Level	Show-dialog	Show status	ignore
err_29035	board_error	red	true	true	false
err_29036	board_error	red	true	true	false
err_29037	board_error	red	true	true	false
err_19038	board_error	yellow	false	true	false
err_19039	board_error	yellow	false	true	false
err_19040	board_error	yellow	false	true	false
err_19041	board_error	yellow	false	true	false
err_19042	board_error	yellow	false	true	false
err_19043	board_error	yellow	false	true	false
err_19044	board_error	yellow	false	true	false
err_19045	board_error	yellow	false	true	false
err_19046	board_error	yellow	false	true	false
err_19047	board_error	yellow	false	true	false
err_19048	board_error	yellow	false	true	false
err_19049	board_error	yellow	false	true	false
err_19050	board_error	yellow	false	true	false
err_19051	board_error	yellow	false	true	false
err_19052	board_error	yellow	false	true	false
err_29053	board_error	red	true	true	false
err_29054	board_error	red	true	true	false
err_29055	board_error	red	true	true	false
err_29056	board_error	red	true	true	false
err_29057	board_error	red	true	true	false
err_29058	board_error	red	true	true	false
err_29059	board_error	red	true	true	false
err_29060	board_error	red	true	true	false
err_29061	board_error	red	true	true	false
err_29062	board_error	red	true	true	false

Error code	Description	Comment
err_10001		
err_10002		
err_10003		
err_10004		
err_10005		
err_10006	Max. OPR HRS total reached	no comment
err_10007	Max. OPR HRS UV reached	no comment
err_10008	Max. OPR HRS el.ballast reached	no comment

Error code	Description	Comment
err_10009	Max. OPR HRS valve 1 reached	no comment
err_10010	Max. OPR HRS valve 2 reached	no comment
err_10011	Max. OPR HRS valve 3 reached	no comment
err_10012	Max. OPR HRS p. tube 1 reached	no comment
err_10013	Max. OPR HRS p. tube 2 reached	no comment
err_10014	Max. OPR HRS pump 1 reached	no comment
err_10015	Max. OPR HRS pump 2 reached	no comment
err_10016	Max. OPR HRS total soon	no comment
err_10017	Max. OPR HRS UV soon	no comment
err_10018	Max. OPR HRS el.ballast soon	no comment
err_10019	Max. OPR HRS valve 1 soon	no comment
err_10020	Max. OPR HRS valve 2 soon	no comment
err_10021	Max. OPR HRS valve 3 soon	no comment
err_10022	Max. OPR HRS pump tube 1 soon	no comment
err_10023	Max. OPR HRS pump tube 2 soon	no comment
err_10024	Max. OPR HRS pump 1 soon	no comment
err_10025	Max. OPR HRS pump 2 soon	no comment
err_10026	Closed Input	no comment
err_10027	Cond.1 wrong Voltage	no comment
err_10028	Cond.1 wrong Frequence	no comment
err_10029	Cond.1 wrong Cell.const	no comment
err_10030	Cond.1 wrong Temp.coeff	no comment
err_10031	Cond.1 wrong Setup time	no comment
err_10032	Cond.1 wrong Hold.time	no comment
err_10033	Cond.2 wrong Voltage	no comment
err_10034	Cond.2 wrong Frequence	no comment
err_10035	Cond.2 wrong Cell.const	no comment
err_10036	Cond.2 wrong Temp.coeff	no comment
err_10037	Cond.2 wrong Setup time	no comment
err_10038	Cond.2 wrong Hold.time	no comment
err_10039	Max. TOC value exceeded	no comment
err_10040	Max. Temp 1 value exceeded	no comment
err_10041	Max. Temp 2 value exceeded	no comment
err_10042	Max. Cond 1 value exceeded	no comment
err_10043	Max. Cond 2 value exceeded	no comment
err_10044	Max. pressure value exceeded	no comment
err_10045	Max. Sys Temp value exceeded	no comment

Error code	Description	Comment
err_10046	Warnlevel TOC Port 1 exceeded	no comment
err_10047	Warnlevel TOC Port 3 exceeded	no comment
err_10048	Warnlevel TOC Port 4 exceeded	no comment
err_10049	Warnlevel Temp 1 exceeded	no comment
err_10050	Warnlevel Temp 2 exceeded	no comment
err_10051	Warnlevel Cond 1 exceeded	no comment
err_10052	Warnlevel Cond 2 exceeded	no comment
err_10053	Warnlevel pressure exceeded	no comment
err_10054	Warnlevel Sys Temp exceeded	no comment
err_10055	Alarmlevel TOC Port 1 exceeded	no comment
err_10056	Alarmlevel TOC Port 3 exceeded	no comment
err_10057	Alarmlevel TOC Port 4 exceeded	no comment
err_10058	Alarmlevel Temp 1 exceeded	no comment
err_10059	Alarmlevel Temp 2 exceeded	no comment
err_10060	Alarmlevel Cond 1 exceeded	no comment
err_10061	Alarmlevel Cond 2 exceeded	no comment
err_10062	Alarmlevel pressure exceeded	no comment
err_10063	Alarmlevel Sys Temp exceeded	no comment
err_10064	RTC time deviates from system time	no comment
err_10065	One or more values are simulated	Values are set to simulate
err_20001	File rights corrupt	Check the rights file hash
err_20002	File users corrupt	Check the users file hash
err_20003	File settings corrupt	Check the settings file hash
err_20004	State not found	State idle check loop
err_20005	Connection lost	Connection to IO Board failed
err_20006	Fluidic leakage	Leakage check
err_20007	UV radiation low	UV Check
err_20008	No Hardware Connection	Connection to IO Board on start
err_20009	String Input Error	
err_20010	Empty string	
err_20011	Not displayable ASCII	
err_19001	Incorrect RS232 address	
err_19002	Probably missing "ROGER"	
err_19003	Error in the request telegram	
err_19004	No command recognised	
err_19005	String is not complete for decoding	
err_29006	Timeout while writing to I2C (RX)	

Error code	Description	Comment
err_29007	Timeout while writing to I2C (TX)	
err_29008	Timeout communication	
err_29009	PC has sent incorrect parameters	
err_19010	The addressed port is reserved	
err_29011	Device responds with incorrect parameter	
err_29012	DA Error	
err_29013	Error initialising DA1	
err_29014	Error initialising DA2	
err_29015	Error initialising DA3	
err_19016	Incorrect gain value (DA)	
err_19017	Incorrect channel (DA)	
err_19018	Transmission error (DA)	
err_19019	Out of range (DA), try to set value < 0	
err_29020	Init IC1 failed	
err_29021	Init IC2 failed	
err_29022	Error when setting a byte of DigOut1	
err_29023	Error when setting a byte of DigOut2	
err_29024	Error when setting a port of DigOut1	
err_29025	Error when setting a port of DigOut2	
err_19026	Port out of range, only 1-16 possible	
err_19027	Only IC1 or IC2 possible	
err_19028	Value outside the value range (DO)	
err_19029	Port out of range, only 1-16 possible	
err_29030	Init IC3 failed	
err_29031	Error when initialising AD1	
err_29032	Error when initialising AD2	
err_29033	Error when initialising AD3	
err_29034	Error when initialising AD4	
err_29035	Sampling value not yet ready	
err_29036	Sampling value request error	
err_29037	Sampling value set error	
err_19038	Incorrect AD-Channel	

Error code	Description	Comment
err_19039	Incorrect sampling rate	
err_19040	Incorrect amplification	
err_19041	Incorrect damping	
err_19042	Incorrect duty value transferred	
err_19043	Incorrect I-Out channel	
err_19044	I-Value out of range	
err_19045	Frequency out of range	
err_19046	Incorrect addressing in the EEPROM	
err_19047	Incorrect number of bytes to be saved	
err_19048	Incorrect motor selection	
err_19049	Incorrect choice of direction for motor	
err_19050	Incorrect speed for motor	
err_19051	Automatic stop too large	
err_19052	Function has been disabled	
err_29053	TimeOut on reception of RS485	
err_29054	Error setting date/time	
err_29055	Init IC3 failed	
err_29056	Error when setting a byte of DigOut3	
err_29057	Error when sending to WD	
err_29058	Error when querying WD	
err_29059	Error when sending to RT	
err_29060	Error when querying RT	
err_29061	Error when sending to TA	
err_29062	Error when querying TA	

Chapter 13

Maintenance

13.1. Maintenance schedule

Maintenance intervals depend heavily on the usage behavior of the appliance. The intervals must therefore be adapted to the specific requirements, but care must be taken to carry out this maintenance work regularly.

Table 13.1.1.: Maintenance plan for miniTOC

Interval	Maintenance work
Before each calibration	Replacing the calibration solution
every 6 months	Replacing the pump hose in case of device with single pump
every 12 months	Replacing the pump hose in case of device with two pumps
every 12 months	Replacing the UV-lamp
every 24-36 months	Replacing the UV reactor Only to be carried out by service personnel authorized by membraPure GmbH.
every 24-36 months	Replacing the UV-lamp ballast Only to be carried out by service personnel authorized by membraPure GmbH.
every 36-48 months	Replace pump head Only to be carried out by service personnel authorized by membraPure GmbH.

13.2. Maintenance work

Improper connection can lead to injury or death! Before starting maintenance work, ensure that there is no voltage on any cable. Disconnect the appliance from the mains and pull out the mains plug.

13.2.1. Peristaltic pump

1. Stop measurements
2. Shut down system
3. Switch off the appliance with the mains switch by setting the cam switch to the "0" position
4. Unplug the appliance from the mains.
5. Open the appliance with the enclosed key.
6. Open pump head
7. Remove the hose and open the hose clamps.
8. Pull new hose over the sample lines
9. Secure both ends with hose clamps
10. Insert the pump hose into the pump head
11. Close pump head
12. Reset operating hours (see chapter 10.2.2).

13.2.2. UV lamp and reactor

Improper handling may cause injury to the skin and eyes. Improper handling can lead to skin and eye injuries.

- Switch off the device before working on the reactor (see chapter 9.12) and disconnect from the power supply.
- The UVC lamp and reactor are separate units; replacing the UV lamp does not involve replacing the reactor!
- Damaged reactors must be replaced
- Never operate the reactor in a shielded state
- Dispose of broken or defective UV lamps as hazardous waste as they contain mercury.
- After switching off the appliance, open the door and wait 30 minutes.

Changing the UV lamp

1. Switch off the system and disconnect from the mains
2. Open the front door using the switch cabinet key
3. Open the cover by turning the four knurled nuts anticlockwise
4. Remove the UV cover.
5. Disconnect the plug connections for the electrical connection of the UV lamp and remove the supply cables from the clamps provided for this purpose
6. Loosen the retaining screw of the fluidic plate. Fold the entire fluidic plate slightly forwards
7. Pull the UV lamp out of the glass spiral reactor upwards - **Handle with care - risk of breakage!**
8. Push the new UV lamp into the center tube of the glass coil reactor until the UV lamp reaches the lower stop. When inserting the UV lamp, the lower connecting cables and the plug connector must first be threaded through the glass coil reactor. - **Handle with care - risk of breakage!**
9. Return the front panel of the appliance to its original position and secure it with the corresponding screws
10. Plug the connections of the UV lamp into the corresponding plug. Attach the UV lamp supply cable to the retaining clips provided.
11. Place the UV cover on the corresponding bolts. Ensure that the sample tubes pass through the recess on the left-hand side.
12. Fasten the UV cover by turning the four knurled screws clockwise.
13. Then reset the operating hours of the UV lamp (see chapter 10.2.1)

13.2.3. Changing the fluidic tubing

To ensure that the hose connectors fit correctly, we recommend cutting the hoses to length using a special hose cutter to ensure a straight cut surface and avoid deformation of the hose ends.

This chapter describes the replacement / repair of fluidic hoses on the device. The hoses used are made of FEP (fluoroethylene propylene) and are resistant to aqueous media in the pure and ultrapure water range.

UNF fittings made of PEEK (polyetheretherketone) are used for the connection. The connection consists of the fitting and a ferrule. The ferrule is mechanically deformed during each installation and must be replaced each time a repair is carried out. If this is not done, it cannot be guaranteed that the device will function properly.


Figure 13.2.1.: Ferrule [1], fitting [2], hose [3]

1. Loosen the fitting of the hose to be replaced by turning it counterclockwise
2. Remove the hose section from the housing
3. Determine the total length of the hose
 - The length is measured from ferrule to ferrule
4. Cut the new hose to the measured length using a suitable hose cutter
 - Make sure the cutting edge is straight!
5. First push the fitting onto the hose


Figure 13.2.2.: Assembly of the fluidic connection

6. Then place the ferrule on the hose. Pay attention to the orientation!
7. Then screw the new hose into the UNF connection.
8. Tighten the renewed connection hand-tight.


Figure 13.2.3.: Installing the connection (1)


Figure 13.2.4.: Installing the connection (2)

13.3. Decommissioning

Risk of injury and infection through contact with the medium!

- Before loosening the hoses, make sure that no action, e.g. pumping a sample, is running or about to start.
- Absorb any leaking reagents with a disposable cloth and rinse with clean water. Dry cleaned areas with a cloth.

1. Navigate to the main menu and click on Direct Direct
2. Remove all connected solutions from the ports (see chapter 4.2).
3. Connect a bottle of pure or ultrapure water.
4. Flush the system for 30 min.
5. Rinse all reagents and standard channels
6. Remove the bottle with pure or ultrapure water
7. Switch on the pump to empty the sample channels.
 - Open each valve for approx. 30s.
8. Close the ports with connection caps

Chapter 14

Accessories

Chapter 15

Certificates

Chapter 16

FAQ

Chapter 17

Technical data

17.1. Input

Measured variable	TOC
Measuring range	0,3 - 1500 µg/L (ppb)
Input signal	trigger control input (option)

17.2. Output

Output signal	TOC, input conductivity
Failure signal	relais signal (Alarms & TOC limit)
Limit value	Max. 500 Ω [Ohm]
Basic system	2 × 4..20 mA Output according to NAMUR NE43 (TOC & Cond.)
Pro system	up to 4 × 4..20 mA Output according to NAMUR NE43 (TOC & Cond.)

17.3. Performance

Measuring range	TOC
Detection crane (LOD)	0,1 µg/L (ppb)
Response time t_{90}	180s
Number of measuring channels	1 online + 1 for calibration or SST (BASIS); 3 for calibration or SST, multiline up to 3 loops (PRO)
Sample requirement	Standard: 7 mL/min
UV reactor	Hg-Quarz-UV-lamp, range of wavelength 185-254 nm
Maintenance interval	see chapter 13.

17.4. Surroundings

Ambient temperature	10°C - 45°C
Storage temperature	2 – 55 °C (35 – 131 °F)
Relative humidity	10 % - 95 %non-condensing
Protection class	IP 54(standard device)
EMC	Emitted interference and interference immunity according to EN 61326- 1:2013 Class A for industrial areas
Electrical safety	2
Degree of soiling	2

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